## INSTRUCTIONS FOR USE

### INTENDED USE

CD-Chex CD103 Plus® is intended to be used as quality control material for evaluating intracellular and surface antigens CD103, CD30, CD38, CD56, CD138, and Cytoplasmic Lambda, with monoclonal antibody binding by flow cytometry. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry, they provide reference values for abnormal peripheral blood leukocytes. CD-Chex CD103 Plus is designed for use on BD Biosciences and Beckman Coulter® flow cytometry systems. **This product and the markers provided on the assay have not been cleared by the U.S. Food and Drug Administration for In Vitro Diagnostic use. This product and the values provided are For Research Use Only. Not for use in diagnostic procedures.**

### SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry is a complex process. CD-Chex CD103 Plus is designed to represent abnormal peripheral blood leukocytes that possess surface antigens and intracellular antigens detectable with fluorescent monoclonal antibodies by flow cytometry<sup>1,2,3</sup>. When processed in the presence of CD45, two distinct populations are clearly visible. Stable antigens on the CD45+ abnormal lymphocyte population include CD103 and CD38. Conversely, an abnormal CD45- cell population present in CD-Chex CD103 Plus provides surface CD30, CD38, CD56, CD138 and intracellular Lambda. CD-Chex CD103 Plus is a positive procedural assayed control used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance to provide consistent and reliable quality control measurements<sup>4</sup>.

### REAGENTS

CD-Chex CD103 Plus contains stabilized human blood and cells of human origin in a preservative medium.

### PRECAUTIONS

1. CD-Chex CD103 Plus is For Research Use Only. Not for use in diagnostic procedures.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at [www.streck.com](http://www.streck.com) for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates the use of the product.
5. CD-Chex CD103 Plus should not be used as a calibrator.
6. SDS can be obtained at [www.streck.com](http://www.streck.com), by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

CD-Chex CD103 Plus is stable through the expiration date when stored at 2 °C to 10 °C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. DO NOT FREEZE.


### INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:


1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of the CD-Chex CD103 Plus. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

### INSTRUCTIONS FOR USE

1. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. Mixing Procedure (**mechanical mixing by vortex or rotator is not recommended**):  
For a video demonstration, visit [www.streck.com/mixing](http://www.streck.com/mixing).
  - a. Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- b. Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.


- c. Aliquot immediately after mixing.
- d. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.

Note: Vials stored for an extended period of time may require extra mixing.

4. Return control reagent to refrigeration after sampling to ensure maximum open-vial stability.
  5. Add CD45 and other monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.
- Note: A negative staining control is recommended due to heterogeneous/dim expression of selected assayed parameters.
6. Incubate according to antibody manufacturer's instructions.
  7. Add recommended amount of RBC lysing agent and follow manufacturer's instructions.
  8. Analyze by flow cytometry following the gating strategy outlined in the Gating section.

### GATING

The most common gating strategy used in neoplastic cell assessment is to gate on the abnormal cells and then determine the CD marker percent positivity using a negative staining control<sup>1-4</sup>. Abnormal cells are generally located on a FSC/SSC plot or a CD45/SSC plot<sup>1-4</sup>. CD45 must be used to separate two cell populations. Debris should be excluded from the cell gate to obtain percent recovery values within assay range.

**FIGURE 1. Gating methods used to identify two populations in CD-Chex CD103 Plus.**

In flow cytometry, abnormal cells can routinely be identified using a FSC/SSC plot (A. Rectangle Gate). However, this gating strategy does not sufficiently allow for the detection of two populations. **The gating strategy employing CD45/SSC (B. Oval Gates) is required for optimal recovery of assayed parameters.** This method easily separates the CD45+ and CD45- cell populations present in CD-Chex CD103 Plus.

### INTRACELLULAR MARKER INSTRUCTIONS FOR USE

Follow steps 1-3 above. CD-Chex CD103 Plus is stabilized cellular material. Therefore, the fixation step (Reagent 1 or Reagent A) used prior to permeabilization in commercially available intracellular staining kits is not required. Use of the fixation reagent will result in sub-optimal recovery.

### LIMITATIONS

To differentiate between the two cell populations present in the CD-Chex CD103 Plus, CD45 must be used during staining. Failure to follow these instructions will result in values outside of the established assay ranges (see gating section).

### EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on properly compensated flow cytometers<sup>4</sup>. The assay values are obtained using common flow cytometry reagents. See the assay for additional limitations or specific instructions for reagents.

**Use of alternate gating strategies not specified in these instructions may result in values outside the published range.** Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different reagents, laboratory protocols, instrument calibration, maintenance, and operator technique. Data collected from interlaboratory quality control programs can be used as a cumulative approach when calculating ranges.

### REFERENCES

1. Gorczyca, W. Flow cytometry in neoplastic hematology, morphologic-immunophenotypic correlation - Second Edition. Informa Healthcare. 2010.
2. Cheria, S., Wood, B. Flow cytometry in evaluation of hematopoietic neoplasms; a case-based approach. 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation - Second Edition. 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematology cells; Approved Guideline - Second Edition. 2009.

### QUALITY CONTROL PROGRAM

Streck offers **STATS®**, an interlaboratory quality control program, to all customers at no charge. For more information, contact the **STATS** Department at 800-898-9563 or [statsdata@streck.com](mailto:statsdata@streck.com). Additional information can be found at [www.streck.com](http://www.streck.com).

### ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [www.streck.com](http://www.streck.com).

See [www.streck.com/patents](http://www.streck.com/patents) for patents that may be applicable to this product.

The brand and product names of the instruments are trademarks of their respective holders.

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## RÉFÉRENCES

- Gorczyca, W. Flow cytometry in neoplastic hematology, morphologic-immunophenotypic correlation - Second Edition. Informa Healthcare. 2010.
- Cheria, S., Wood, B. Flow cytometry in evaluation of hematopoietic neoplasms; a case-based approach. 2012.
- Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation - Second Edition. 2007.
- Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematology cells; Approved Guideline - Second Edition. 2009.

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## QUELLANGABEN

1. Gorczyca, W. Flow cytometry in neoplastic hematology, morphologic-immunophenotypic correlation - Second Edition. Informa Healthcare. 2010.
2. Cheria, S., Wood, B. Flow cytometry in evaluation of hematopoietic neoplasms; a case-based approach. 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation - Second Edition. 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematomatoid cells; Approved Guideline - Second Edition. 2009.

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1. Gorczyca, W. Flow cytometry in neoplastic hematology, morphologic-immunophenotypic correlation - Second Edition. Informa Healthcare. 2010.
2. Cheria, S., Wood, B. Flow cytometry in evaluation of hematopoietic neoplasms; a case-based approach. 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation - Second Edition. 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematomalymphoid cells; Approved Guideline - Second Edition. 2009.

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## REFERENSER

- Gorczyca, W. Flow cytometry in neoplastic hematology, morphologic-immunophenotypic correlation - Second Edition. Informa Healthcare. 2010.
- Cheria, S., Wood, B. Flow cytometry in evaluation of hematopoietic neoplasms; a case-based approach. 2012.
- Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation - Second Edition. 2007.
- Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematology lymphoid cells; Approved Guideline - Second Edition. 2009.

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