INSTRUCTIONS FOR USE

INTENDED USE

Cell-Chex® Auto is an assayed whole blood control for evaluating the accuracy and precision of hematology instruments that measure blood cell counts in patient body fluid samples.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, Cell-Chex Auto will provide values within the expected range indicated on the assay sheet.

REAGENTS

This product contains stabilized human red blood cells and white blood cells in a preservative medium.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
- Control products are not to be used as calibrators.
- SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

Cell-Chex Auto is stable through the expiration date when stored at 2 °C-10 °C. After opening, Cell-Chex Auto is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C-10 °C.

INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:

- Review control product package insert and the operating procedure of the instrument.
- Check expiration date of the product on the vial. Discard out-dated products.
- Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

- Remove a vial of the control from the refrigerator and warm to room temperature (18 °C-30 °C) for 15 minutes before use.
- To mix: **(Do NOT mix mechanically or vortex.)**
For a video demonstration, visit streck.com/mixing.
a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
 - Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).
- Begin with a system rinse on all systems before running Cell-Chex Auto to reduce carryover.
NOTE: It is critical for all systems that the background be low prior to running body fluid controls. Carryover from previous samples is significant compared to the cell counts contained in a body fluid sample and could result in erratically high cell counts.

SPECIAL INSTRUCTIONS FOR SELECTED INSTRUMENTS

- Abbott® CELL-DYN Ruby®**
 - Begin with a system rinse.
 - Enter Sample I.D. Select "Patient" as the specimen type and "CBC & RRBC" as the test selection.
 - WOC and NOC values for WBC are provided on the assay. Report desired WBC value(s) per laboratory body fluid protocol.
 - Start with the lowest level of control during each control analysis period to reduce carryover.
- Beckman Coulter® UniCel® DxH® 600/800**

A. QC Setup using 2D Barcode

 - Select: **Menu/Setup/Quality Control**
 - Select: **New Control from Barcode**
 - Scan 2D Barcode
 - Verify files have uploaded correctly and select "OK"

B. Running Cell-Chex Auto

- Select: **Single-tube presentation**
 - Scan vial barcode with handheld barcode scanner**
 - Verify DxH control lot number code is displayed in the "Specimen Identifier" box
 - Check the "DxH Control" box
 - Click the cursor inside the "Specimen Identifier" box and press Enter
 - Click "OK" on the DxH "Test Order Recognition" box
 - Place tube into manual aspiration holder
- Siemens Healthcare Diagnostics ADVIA 120/2120/2120i**
 - Begin with a system rinse.
 - Start with the lowest level of control during each control analysis period to reduce carryover.
 - Run manually in the open mode.
 - WBC-P values are provided on the assay.
 - Sysmex XE-2100, XT-1800i, XT-2000i**
NOTE: Control files do not support Diff Channel WBC values. To locate the WBC#(DIFFch) parameter: from the sample run, go to the service tab and select the diff information in the drop-down bar. The WBC#(DIFFch) is toward the middle of the screen.
 - Begin with a system rinse.
 - Start with the lowest level of control during each control analysis period to reduce carryover.
 - Baso Channel assay value (reported as WBC x 10⁹/L) and Diff Channel values for WBC are provided on the assay. Report desired WBC value(s) per laboratory body fluid protocol.
 - Cell-Chex Auto may be run in open or closed mode.
 - Sysmex XE-5000, XT-4000i** - Run Cell-Chex Auto in the Body Fluid Mode. Follow the instrument manual instructions for running samples in the Body Fluid Mode.

After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both vial and cap before replacing cap and returning to refrigeration.

LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with this product. The white blood cell components do not simulate white blood cell morphology.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

QUALITY CONTROL PROGRAM

Streck offers **STATS®**, an interlaboratory quality control program, to all customers at no charge. For more information, contact the **STATS** Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

All product names, logos, brands, and marks are property of their respective owners.

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