## INSTRUCTIONS FOR USE

### INTENDED USE

Cal-Chex® is a quantitative calibrator intended to be used for calibrating automated multi-parameter hematology instruments for Red Blood Cell Count (RBC), White Blood Cell Count (WBC), Hemoglobin (HGB), Mean Cell Volume (MCV), Hematocrit (HCT), Mean Platelet Volume (MPV), and Platelet Count (PLT).

### SUMMARY AND PRINCIPLES

Hematology analyzers require regular calibration in order to produce accurate results on patient samples. Calibration can be accomplished by transferring information to the analyzer through fresh blood samples which have been assayed by reference methods. A more direct and convenient approach is to use a calibrator material with System Specific Values (SSV) assigned.

Streck Calibrators are stable suspensions of red blood cells, white blood cells and platelets. Assigned values are derived from replicate analysis on whole-blood calibrated hematology analyzers. Users analyze the calibrator on their instruments and compute calibration factors by comparing recovered values and assigned values. These factors provide the basis for making adjustments to the instrument.

### REAGENTS

This whole blood reagent may contain any or all of the following: stabilized human or animal red blood cells, human, animal or simulated white blood cells and a platelet component in a preservative medium.

### PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at [streck.com](http://streck.com) for specific FDA required blood tests.
3. This product should not be disposed of in general waste, but should be disposed of with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. SDS can be obtained at [streck.com](http://streck.com), by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY


Cal-Chex is stable through the expiration date when stored at 2 °C to 10 °C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C.

### INDICATION OF PRODUCT DETERIORATION

Gross hemolysis (darkly colored supernatant) may be indicative of product deterioration or damage. If this is apparent, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

### INSTRUCTIONS FOR USE

1. Perform instrument start-up and routine cleaning procedures as defined in instrument manufacturer's operator's manual (i.e., blood sampling valve, counting apertures).
2. Remove vials of calibrator from refrigeration and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. To mix: **(Do NOT mix mechanically or vortex.)**  
For a video demonstration, visit [streck.com/mixing](http://streck.com/mixing).
  - a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
  - d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).
4. Refer to the instrument manufacturer's operator manual for detailed instructions on priming the instrument with normal patient samples (if applicable), verifying instrument repeatability according to their specifications, executing the instrument auto-calibration process, and performance of calibration runs.
  5. Calculate the mean value for each parameter.
  6. Compare the results of your calculations to the values listed on the assay for your instrument.
    - a. If the difference between your recovered mean values and the system specific values are less than the listed tolerance limits, the instrument does not require calibration.
    - b. If the difference is greater, calibrate using the system specific values.
  7. Calibration of the specific parameter(s) should be done in accordance with the procedure in your instrument manual. Calibration may not be required for all parameters.
  8. Refer to instrument manufacturer's operator manual to determine if additional testing is required to verify calibration changes (if applicable) and/or verify calibration.
  9. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both the vial and cap before replacing cap and returning to refrigeration.
  10. For further assistance, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

### EXPECTED RESULTS

Assay values are derived from replicate testing on instruments operated and maintained according to manufacturer's instructions.

The following reference methods are used to whole blood calibrate instruments prior to system specific value assignment as indicated in CLSI H26-A2.<sup>1</sup>

- WBC** - A series of 1:500 dilutions are made using class A glassware. The lytic agent is placed in the initial dilution flask before diluting to volume. The diluting agent is an isotonic non-interfering solution. The samples are counted on a single aperture volume displacing particle counter.
- RBC** - A series of 1:50,000 dilutions are made using class A glassware. The diluting agent is an isotonic non-interfering solution. The samples are counted on a single aperture volume displacing particle counter.
- Hb** - Hemoglobin concentration is determined by converting hemoglobin to hemoglobincyanide (HiCN) and measuring absorbance at 540nm according to CLSI H15-A3<sup>2</sup> and ICSH recommendations. Hemoglobin concentration is calculated using millimolar absorption of 11.0.
- Hct** - Microhematocrit values are done in replicate on each sample, with capillary tubes filled and centrifuged according to the CLSI H07-A3<sup>3</sup> document. K<sub>2</sub>EDTA is used as the anticoagulant for the collection of fresh specimens. The packed cell volume (hematocrit) is read directly using a precision metric scale. No correction is made for trapped plasma.
- Plt** - Platelet counting by the RBC/platelet ratio method (ICSH) is utilized as the reference method.

### REFERENCES

1. Clinical Laboratory Standards Institute, H26-A2. Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
2. Clinical Laboratory Standards Institute, H15-A3. Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard - Third Edition.
3. Clinical Laboratory Standards Institute, H07-A3. Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

### ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at [streck.com](http://streck.com).

### Rx Only

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [streck.com](http://streck.com).

See [streck.com/patents](http://streck.com/patents) for patents that may be applicable to this product.

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**RÉFÉRENCES**

1. Clinical Laboratory Standards Institute, H26-A2. Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
2. Clinical Laboratory Standards Institute, H15-A3, Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard - Third Edition.
3. Clinical Laboratory Standards Institute, H07-A3, Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

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## QUELLENANGABEN

- Clinical Laboratory Standards Institute, H26-A2. Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
- Clinical Laboratory Standards Institute, H15-A3, Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard - Third Edition.
- Clinical Laboratory Standards Institute, H07-A3, Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

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