INSTRUCTIONS FOR USE

INTENDED USE

A1c-Cellular[®] is a bi-level, whole blood based, assayed control material for evaluating the accuracy and precision of HbA1c procedures.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated and semi-automated procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of the instrument. A1c-Cellular's intact RBC formulation allows users to verify all steps of analysis procedures, including lysing of the RBC. A1c-Cellular is for use by laboratory professionals.

REAGENTS

A1c-Cellular contains stabilized human red blood cells in a preservative medium. No additional materials are provided or required.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product, unless pretreatment is recommended by individual instrument manufacturers.
- Control products are not to be used as calibrators.
- SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.
- This product contains a substance toxic to reproduction, at non-hazardous levels, and a preservative which may result in a sensitization or an allergic reaction. All users are to wear standard laboratory PPE.

STORAGE AND STABILITY

A1c-Cellular is stable through the expiration date when stored unopened at 2 °C to 10 °C. After initial opening, A1c-Cellular is stable for 30 days when stored at 2 °C to 10 °C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration; however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:


- Review control product package insert and the operating procedure of the instruments.
- Check expiration date of the product on the vial. Discard outdated products.
- Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

- A1c-Cellular should only be used with the instruments listed on the assay.
- Remove vials of A1c-Cellular from the refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
- To mix: **(Do NOT mix mechanically or vortex.)**

For a video demonstration, visit streck.com/mixing.

- Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
 - Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).
- Refer to the instrument manual for the system in use for analyzing control materials.
 - After sampling, return to refrigeration for maximum open-vial stability.

Barcodes on the cap-pierceable vials are intended for use with Bio-Rad[®] VARIANT II TURBO and Tosoh[®] G8. To use the barcode, align the barcode on the vial(s) with the barcode reader on the instrument. The barcode will label the control with "LC" for Level 1 and "HC" for Level 2. The Bio-Rad Variant instrument can use these barcodes to transfer the data into a QC log file. In order to not use the barcode function, face the barcode away from the instrument's barcode reader. Report any serious incident to Streck and appropriate regulatory entities including the competent authority of the member state in which the user and/or patient is established as applicable.

LIMITATIONS

Do not use the product past the expiration date or if there is evidence of gross hemolysis or microbial contamination.

EXPECTED RESULTS AND THEIR DERIVATION

The mean assay values provided are derived from replicate analysis. The values are lot specific and were performed by the control manufacturer and/or independent laboratories using the same lot. The mean

values recovered on your system should be within the stated expected range. Good laboratory practice includes establishment of a laboratory mean and standard deviation by replicate analyses. Assay ranges are provided only as a guide.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2¹ it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

REFERENCES

- Clinical and Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers; Approved Standard - Second Edition.

QUALITY CONTROL PROGRAM

Streck offers STATS[®], an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-843-0912 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

CHANGES FROM PREVIOUS VERSION

Updated to meet requirements of Regulation (EU) 2017/746.

EU SUMMARY OF SAFETY AND PERFORMANCE (SSP)

SSP is available on the EUDAMED website (ec.europa.eu/tools/eudamed/).

All product names, logos, brands, and marks are property of their respective owners.

See streck.com/patents for patents that may be applicable to this product.

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