## MANUAL VIRAL EXTRACTION KIT

### INSTRUCTIONS FOR USE

#### INTENDED USE

The Streck® Viral Extraction Kit is intended for use in isolating intact virus particles and associated viral nucleic acids (RNA) contained in clinical samples for use in Molecular Diagnostic Assays. Utilizing Ceres Nanotrap® Microbiome A Particles and Streck developed buffers, viral RNA extraction and purification is completed in less than 30 minutes. Eluted samples can then be utilized for downstream analysis or stored for later use.

#### SUMMARY AND PRINCIPLES

Nanotrap Microbiome A Particles isolate and concentrate intact viruses and virus-like particles (VLPs) from samples in an affinity matrix. The Nanotrap Particles, attached viruses and VLPs are bound by a strong magnetic attraction to the magnets such as the DynaMag™-2 rack. The viruses are lysed, and the nucleic acids bound to the Nanotrap. Nucleic acid samples are then washed and eluted into a sample vial for use in clinical or research assays.

#### PRECAUTIONS

- For In Vitro Diagnostic Use.
- Nanotrap Microbiome A Particles are intended for use as supplied. Do not dilute or add other components to the Nanotrap Microbiome A Particles.
- SDS can be obtained at [streck.com](http://streck.com), by calling 800-843-0912 or by calling your local supplier.
- Dispose of used reagents, viral debris and consumables as hazardous waste according to established laboratory procedures.
- Always use DNase/RNase free plasticware/reagents and aerosol-barrier pipet tips.

#### STORAGE AND STABILITY

- Store between 2 °C - 25 °C. DO NOT FREEZE.
- When stored at the proper temperature range, kit contents are stable through the expiration date.

#### INDICATIONS OF PRODUCT DETERIORATION

- If Nanotrap Microbiome A Particles have previously been frozen, large visible aggregates may form once the material is thawed. If aggregates are present, product should not be used.
- If frozen, do not use.

#### MATERIALS PROVIDED IN THE STRECK VIRAL EXTRACTION KIT

- Buffer A – Lysis/Binding Solution
- Buffer B – Wash/Elution Buffer
- Nanotrap Microbiome A Particles

#### MATERIALS REQUIRED BUT NOT SUPPLIED

- PCR clean microfuge tubes with locking lids
- PCR clean, filter tips for pipets
- Adjustable Volume Pipets
- DynaMag-2 Magnet (ThermoFisher Scientific Catalog #12321D) (or equivalent)
- Heat block capable of holding 1.5 – 2.0ml microfuge tubes
- Freezer capable of storing samples at -80 °C (If not analyzing the purified viral RNA immediately)
- Vortex
- Laboratory timer

#### PROTOCOL

- Turn heat block on and set temperature to 95 °C.
- Label a microcentrifuge tube for each sample.
- Add the following to each tube:
  - Add 500ul of Sample.
  - Add 400ul of Buffer A.
  - Vortex Nanotrap Microbiome A Particles for 5 seconds at maximum and add 100ul to each tube.
- Vortex samples for 5 seconds.
- Incubate samples for 5 minutes at Room Temperature.
- Place sample tubes on DynaMag-2 Magnet for a minimum of 1.5 minutes.
- Discard supernatant.
- Add 500ul of Buffer A to each sample.
- Vortex samples for 5 seconds.
- Place sample tubes on DynaMag-2 Magnet for a minimum of 1.5 minutes.
- Remove and discard Buffer A wash.
- Add 500ul Buffer B to each sample.
- Vortex sample for 5 seconds.
- Place tubes back on DynaMag-2 Magnet for a minimum of 1.5 minutes.
- Remove and discard Buffer B wash.
- Remove any wash leftover in cap.
- Add 50ul of Buffer B to each sample.
- Vortex sample for 5 seconds.
- For optimal results, place sample tubes on Heat Block and incubate at 95 °C for 5 to 10 minutes.
- Place samples on DynaMag-2 Magnet for 1.5 minutes.
- Transfer supernatant to new labeled nuclease free microcentrifuge tube.
  - Tip magnet to fully separate Nanotrap Microbiome A Particles from sample during pipetting.
- Proceed immediately to downstream analysis or store RNA at -80 °C.

## AUTOMATED VIRAL EXTRACTION

### INSTRUCTIONS FOR USE

#### INTENDED USE

The Streck Viral Extraction Kit is intended for use in isolating intact virus particles and associated viral nucleic acids (RNA) contained in clinical samples for use in Molecular Diagnostic Assays. Utilizing Ceres Nanotrap Microbiome A Particles and Streck developed buffers, viral RNA extraction and purification is completed in less than 30 minutes. Eluted samples can then be utilized for downstream analysis or stored for later use.

#### SUMMARY AND PRINCIPLES

Nanotrap Microbiome A Particles isolate and concentrate intact viruses and virus-like particles (VLPs) from samples in an affinity matrix. The Nanotrap Particles, attached viruses and VLPs are bound by a strong magnetic attraction to the magnets on the KingFisher Apex. The viruses are lysed, and the nucleic

acids bound to the Nanotrap. They are then washed and eluted into a sample vial for use in clinical or research assays.

#### PRECAUTIONS

- For In Vitro Diagnostic Use.
- Nanotrap Microbiome A Particles are intended for use as supplied. Do not dilute or add other components to the Nanotrap Microbiome A Particles.
- SDS can be obtained at [streck.com](http://streck.com), by calling 800-843-0912 or by calling your local supplier.
- Dispose of used reagents, viral debris and consumables as hazardous waste according to established laboratory procedures.
- Always use DNase/RNase free plasticware/reagents and aerosol-barrier pipet tips.

#### STORAGE AND STABILITY

- Store between 2 °C - 25 °C. DO NOT FREEZE.
- When stored at the proper temperature range, kit contents are stable through the expiration date.

#### INDICATIONS OF PRODUCT DETERIORATION

- If Nanotrap Microbiome A Particles have previously been frozen, large visible aggregates may form once the material is thawed. If aggregates are present, product should not be used.
- If frozen, do not use.

#### MATERIALS PROVIDED IN THE STRECK VIRAL EXTRACTION KIT

- Buffer A – Lysis/Binding Solution
- Buffer B – Wash/Elution Buffer
- Nanotrap Microbiome A Particles

#### MATERIALS REQUIRED BUT NOT SUPPLIED

- KingFisher 96 deep-well tip comb
- KingFisher 96 deep-well plate
- KingFisher 96 well standard plate
- KingFisher 96 deep-well magnetic head
- PCR clean, filter tips for pipets
- Adjustable Volume Pipets
- Freezer capable of storing samples at -80 °C (if not analyzing the purified viral RNA immediately)

#### Protocol for using the Streck Viral Extraction Kit on the KingFisher Apex, which can be downloaded at [streck.com/resources](http://streck.com/resources).

- Prepare KingFisher plates following the table below. Vortex Nanotrap Microbiome A Particles for a minimum of 5 seconds at maximum and add 100ul to each well.

	Plate Type	Plate Label	Reagent Added		Position on Apex
1	96 Deep-Well Plate	Sample	Sample	500ul	1
			Buffer A	400ul	
			Nanotrap	100ul	
2	96 Deep-Well Plate	Wash 1	Buffer A	500ul	2
3	96 Deep-Well Plate	Wash 2	Buffer B	500ul	3
4	96 Well Standard Plate	Elution	Buffer B	50ul	4

- Once plates are prepared, the 96 deep-well tip comb must be placed on top of the sample plate.
- Open the extraction protocol program for the Streck Viral Extraction Kit.
  - The program can also be downloaded from [streck.com/resources](http://streck.com/resources).
- The KingFisher will call for each plate to be loaded onto the deck.
  - Plate will be loaded according to the position in the table above.
- Once all plates are loaded, start program.
  - Program will take approximately 22 minutes.
- Once completed, remove final elution plate.
- Proceed immediately to downstream analysis or store RNA at -80 °C.
- Run UV treatment immediately for 1 hour to decontaminate remaining plates as described in the KingFisher procedures.

#### ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at [streck.com](http://streck.com).

#### TECHNICAL SUPPORT

Please contact Streck Technical Services at +1-402-691-7510 (outside the United States and Canada), 800-843-0912 (United States and Canada) or [technicalservices@streck.com](mailto:technicalservices@streck.com). Report any serious incident to Streck, Inc. and appropriate regulatory entities as applicable.

#### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [streck.com](http://streck.com).

All product names, logos, brands, and marks are property of their respective owners.

See [streck.com/patents](http://streck.com/patents) for patents that may be applicable to this product.

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