

ESR-Chex Plus

Instrument Group Name - **Diesse MINI-CUBE**

Manager Acct Number: 20197
Affiliation: A TEST LAB GROUP
ATT: LAB MANAGER

Due Date: 12/5/2020
Lot Expiration: 3/31/2021

Date Prepared: 2/5/2021

Dear Laboratory Manager,

The attached Manager's Report compares the most recent data submitted by each member of your affiliation. Each member who submitted data for this submission is listed below.

Account #	Name	Location	
08588	TEST LAB A	1234 MAIN STREET	ANYTOWN, NE 12345, US
08589	TEST LAB B	200 1ST STREET SW	ANYTOWN, NE 12345, US

If you have any questions about your report, please contact the Streck STATS[®] department at 800-898-9563.

Key to interpreting your report

Lab #	Lab number for each member within the affiliation
Shift #	Indicates the shift for each lab #
Lab N	The number of data points submitted for each parameter by each lab
Assay Mean	The value printed on the product's assay sheet
Lab Mean, SD, CV	Arithmetic calculations for each lab's data
Group Mean, SD, CV	Arithmetic calculations for all of the members of your group
Total Mean, SD, CV	Arithmetic calculations for total # of customers in instrument group
Group SDI	Measure of the difference between each lab's mean and the group mean
Total SDI	Measure of the difference between each lab's mean and the total mean

Reviewed By: _____

Date: _____

ESR-Chex Plus

Instrument Group Name - **Diesse MINI-CUBE**

Manager Acct Number: 20197

Affiliation: A TEST LAB GROUP

ATT: LAB MANAGER

Due Date: 12/5/2020

Lot Expiration: 3/31/2021

Instrument List

Lab #	ID for Report	Serial #	Instrument Name
08588	A	1-0000	Diesse MINI-CUBE
08589	B	2-0000	Diesse MINI-CUBE

ESR-Chex Plus

Instrument Group Name - Diesse MINI-CUBE

Manager Acct Number: 20197	Due Date: 12/5/2020
Affiliation: A TEST LAB GROUP	Lot Expiration: 3/31/2021
ATT: LAB MANAGER	

Current Submission

Parameter: ESR

Lab #	ID	Shift	Lab N	Assay Mean	Lab Mean	Group Mean	Total Mean	Lab SD	Group SD	Total SD	Lab CV	Group CV	Total CV	Group SDI	Total SDI
Level 1 L1-611															
08588	A	1	38	7	2	2	3	1.0	1.2	2.2	50.0	61.8	82.0	0.0	-0.3
08589	B	1	30	7	2	2	3	1.5	1.2	2.2	75.0	61.8	82.0	0.0	-0.3
Level 2 L2-611															
08588	A	1	38	59	59	60	57	5.0	5.2	5.8	8.5	8.6	10.0	-0.3	0.3
08589	B	1	30	59	62	60	57	5.0	5.2	5.8	8.1	8.6	10.0	0.3	0.8

ESR-Chex Plus


Instrument Group Name - Diesse MINI-CUBE

Manager Acct Number: 20197
Affiliation: A TEST LAB GROUP
ATT: LAB MANAGER

Due Date: 12/5/2020
Lot Expiration: 3/31/2021

ESR

Level 1



ESR

Level 2

